

# EICH-COLON®

Medical machines

# COLON-HYDRO



**EN ISO 13485**

**CE 1015**

The Colonic Irrigator offers new options for the treatment of disorders of the colon:

Edition 06/2012-03

- Repolarized water – causes a positive change in the permeability of the cell membrane
- Activated water – promotes vital processes
- Colour therapy – strengthens the individual constitution
- Resonance therapy – collects valuable information on pathological and toxic processes
- Individual patient vibrations – passes on information contained in urine and sputum
- Integrated instillation of oxygen and infusions – attacks the problem directly at the source
- Appealing modern design – pleases the customer's sense of aesthetics
- Wall-mounted and mobile units – offers the therapist greater flexibility

# Advantages of the **Colonic Irrigator** at a Glance

- **Therapeutic options :**
- **tried and testes treatments can be administered without additional time consuming procedures**
- **The Colonic Irrigator comes as a complete instrument unit with an integrated disinfection device**
- **A thermostat-controlled mixing battery ensures a constant water temperature and permits precise temperature adjustment**
- **The instrument fixtures are made of high-quality stainless steel**
- **The unit includes a cold and hot water filter, guaranteeing clean treatment water**
- **All connections on the unit are standardized to ensure quick hook-up**

## List of indications (extract) :

- Constipation
- Diarrhea
- Diverticulitis
- Asthma
- Colitis
- Parasitic infestations
- Inactive colon
- Intestinal toxemias
- Mycoses
- Chronic diseases
- Psychotherapy
- Immunostimulation
- Allergies
- Crohn´s disease

Filtered, **activated water** at a controlled temperature is injected into the colon in the form of a gentle enema. The colon is filled and then emptied in rapid succession; as a result, even obstinate stool deposits are dislodged and flushed out. The colon treatment is a gentle and absolutely hygienic therapy.

With the **color therapy** option, the vibration spectrum of the individual patient can be transferred to the water and conveyed to internal organs. The storage capability of water and the high water content (80%) of the human body are two reasons why this treatment is so successful.

Because of the many different pathological substance and kind of information located in the intestines, it makes sense to apply the **inverted resonance therapy** directly at the site where many pathological processes begin.

To optimize the beneficial effect of this treatment, **endogenous substances** such as the patients own urine and sputum are put to effective use. These substances are instilled in the intestine in **inverted form, i.e. as physiological vibrations.**

The Colonic Irrigator accomplishes more than the traditional **flushing out of toxins**, a procedure which is in itself of undisputed benefit. During colonic irrigator, it simultaneously administers an individual treatment adapted to the requirements of the individual patient **without an extra time consuming procedure.**

Full attention has been paid to hygienic considerations: the Colonic Irrigator Unit comes with a built-in disinfection device and a set of disposable materials.

**Manufacture - Distribution - Service – Disposable Items**

## **EICH - COLON**

- medical machines -  
D - 72793 Pfullingen  
Kunstmühlestr. 12

Fon: +49-7121-72 441 - Fax: +49-7121-790786 - [www.eichcolon.de](http://www.eichcolon.de) - [mail@eichcolon.de](mailto:mail@eichcolon.de)